


## [KO Validated] Anti-Prmt4/Carm1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q86X55

Item Number: CM0002020



### Introduction

KO-validated rabbit monoclonal antibody against PRMT4/CARM1 (Q86X55) for WB, IP, and ELISA. Reacts with human, mouse, and rat samples. Ideal for chromatin biology, transcription regulation, and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] PRMT4/CARM1 Rabbit mAb                                                                                     |
| Remarks/Alias    | PRMT4; M1                                                                                                                 |
| Species          | Human                                                                                                                     |
| Gene ID (Human)  | 10498                                                                                                                     |
| Gene ID          | 10498                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 486-603 of human PRMT4/CARM1 (NP_954592.1). |
| Source           | Rabbit                                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                                     |
| Application      | WB, IP, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | Q86X55                                                                                                                    |
| Protein Weight   | 66kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |