


Ko Validated Anti-Psme2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9UI46

Item Number: CM0016761



Introduction

KO-validated rabbit polyclonal antibody targeting human PSME2 (PA28 β , REG β). Validated for WB and ELISA with cross-reactivity to human and mouse. Essential tool for immunoproteasome and antigen processing research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PSME2 Rabbit pAb
Remarks/Alias	PA28B; REG β ; PA28 β ; E2
Species	Human
GeneID (human)	5721
GeneID	5721
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-239 of human PSME2 (NP_002809.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UL46
Protein Weight	27kDa
Shipping	Ice bag