

[KO Validated] Anti-Shmt2 Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - P34897

Item Number: CM0008460



Introduction

KO validated anti-SHMT2 rabbit monoclonal antibody for WB, IF, IHC, ELISA. Reacts with human, mouse, rat SHMT2. Ideal for mitochondrial one-carbon metabolism and folate pathway research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SHMT2 Rabbit mAb
Remarks/Alias	GLYA; SHMT; mSHMT; NEDCASB; HEL-S-51e
Species	Human
Gene ID (Human)	6472
Gene ID	6472
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 405-504 of human SHMT2 (NP_005403.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34897
Protein Weight	53 kDa/55 kDa/56 kDa
Shipping	Ice bag