


[KO Validated] Anti-Tlk2 Polyclonal Antibody For Wb & Elisa - Q86Ue8

Item Number: CM0028276



Introduction

Rabbit polyclonal antibody against human TLK2, validated for WB and ELISA. Recognizes human and mouse TLK2, a serine/threonine kinase involved in chromatin assembly and DNA damage response. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TLK2 Rabbit pAb
Remarks/Alias	HsHPK, MRD57, PKU-ALPHA, K2
Species	Human
Gene ID (Human)	11011
Gene ID	11011
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human TLK2 (NP_006843.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86UE8
Protein Weight	88kDa
Shipping	Ice bag