


[Ko Validated] NdrG1 Rabbit Pab For Wb, If/Icc, Elisa - Q92597

Item Number: CM0018779



Introduction

KO validated NDRG1 rabbit polyclonal antibody for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Involved in stress response, tumor suppression, and nerve myelin maintenance. Ideal for cancer and neuropathy research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NDRG1 Rabbit pAb
Remarks/Alias	GC4, RTP, DRG1, NDR1, NMSL, TDD5, CAP43, CMT4D, DRG-1, HMSNL, RIT42, TARG1, PROXY1, G1
Species	Human
GeneID (Human)	10397
GeneID	10397
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 175-394 of human NDRG1 (NP_001128714.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92597
Protein Weight	43kDa
Shipping	Ice bag