

Mixl1 Rabbit Pab For Wb And Elisa - Q9H2W2

Item Number: CM0023033

Introduction

Rabbit polyclonal antibody against human MIXL1, validated for Western blot and ELISA. Detects human, mouse, and rat MIXL1; suitable for developmental biology and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	MIXL1 Rabbit pAb
Remarks/Alias	MIX; MIXL; MILD1; MIXL1
Species	Human
GeneID (Human)	83881
GeneID	83881
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 150-232 of human MIXL1 (NP_114150.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2W2
Protein Weight	25kDa
Shipping	Ice bag