


Monomethyl-Histone H2B-K5 Rabbit Polyclonal Antibody (P62807/Q16778) For Wb, If/Icc, Elisa

Item Number: CM0000235



Introduction

MonoMethyl-Histone H2B-K5 Rabbit pAb detects endogenous monomethylated histone H2B at lysine 5 (P62807/Q16778) in human, mouse, and rat. Suitable for WB, IF/ICC, and ELISA, supporting chromatin and epigenetic studies.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H2B-K5 Rabbit pAb
Remarks/Alias	H2B; H2BE; H2BQ; GL105; H2B.1; H2BFQ; H2BGL105; H2B-GL105; HIST2H2BE; MonoMethyl-Histone H2B-K5
Species	Human
GeneID (human)	H2B
GeneID	3017/8349
Immunogen	A synthetic monomethylated peptide around K5 of human H2BK5me1 (NP_003519.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62807/Q16778
Protein Weight	14kDa
Shipping	Ice bag