

Monomethyl-Histone H3-K79 Rabbit Monoclonal Antibody (Q16695/P68431) - Wb, Chip, Ihc-P, Elisa

Item Number: CM0000351



Introduction

Rabbit monoclonal antibody specific for monomethylated lysine 79 of histone H3 (K79me1). Ideal for WB, DB, IHC-P, ELISA, ChIP, and ChIP-seq. Recognizes human, mouse, and rat H3; predicted wide cross-reactivity.

[Learn More](#)

Parameter	Specification
Product Name	MonoMethyl-Histone H3-K79 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; MonoMethyl-Histone H3-K79
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic monomethylated peptide around K79 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag