

Monomethyl-Histone H3-K9 Rabbit Mab | Q16695 / P68431 - Methylation-Specific Antibody For Wb, Chip, Elisa, And More

Item Number: CM0000411



Introduction

Rabbit monoclonal antibody specific for monomethylated histone H3 at lysine 9 (H3K9me1, SWISS Q16695 / P68431). Validated for WB, IF/ICC, IHC-P, DB, ChIP, and ELISA with broad predicted cross-reactivity (human, mouse, rat).

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K9 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic monomethylated peptide corresponding to residues surrounding K9 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, IHC-P, DB, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	16kDa
Shipping	Ice bag