


Monomethyl-Histone H4-K5 Rabbit Pab (P62805) For If/Icc And Elisa

Item Number: CM0000285



Introduction

Rabbit polyclonal antibody targeting monomethylated Histone H4 at lysine 5 (P62805). Validated in IF/ICC and ELISA, with broad species reactivity. Ideal for epigenetic research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H4-K5 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; MonoMethyl-Histone H4-K5
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	Synthetic peptide, a synthetic monomethylated peptide around K5 of human Histone H4 (NP_003529.1)
Source	Rabbit
Category	Methylated Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag