

Monomethyl-Uhrf1-K385 Rabbit Pab (Q96T88) - Epigenetics Antibody For Wb, Db, Elisa

Item Number: CM0000232



Introduction

Rabbit polyclonal antibody specific for monomethylated UHRF1 at K385 (Q96T88). Validated for Western blot, dot blot, and ELISA in human, mouse, and rat samples. Ideal for DNA methylation and chromatin studies.

[Learn More](#)

Parameter	Specification
Product Name	MonoMethyl-UHRF1-K385 Rabbit pAb
Remarks/Alias	Np95; hNP95; ICBP90; RNF106; TDRD22; hUHRF1; huNp95; MonoMethyl-UHRF1-K385
Species	Human
GeneID (human)	29128
GeneID	29128
Immunogen	Synthetic peptide A synthetic monomethylated peptide around K385 of human UHRF1 (NP_001041666.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96T88
Protein Weight	90kDa
Shipping	Ice bag