

Mouse Lins1 Rabbit Polyclonal Antibody For Wb And Elisa - Lins1

Item Number: CM0016913

Introduction

Rabbit polyclonal antibody targeting mouse LINS1, validated for Western blot and ELISA. Detects endogenous LINS1 in mouse and rat samples; predicts bands at 55 kDa and 85 kDa. Suitable for neurological and developmental research.

[Learn More](#)

Parameter	Value
Product Name	LINS1 Rabbit pAb
Remarks/Alias	Lins; Lins2; Wins2; 2700083B01Rik; LINS1
Species	Mouse
Gene ID (Human)	55180
Gene ID	72635
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse LINS1 (NP_690028.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
Protein Weight	55kDa/85kDa
Shipping	Ice bag