


Mrps16 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y3D3

Item Number: CM0030118



Introduction

Rabbit polyclonal antibody against mitochondrial ribosomal protein MRPS16. Validated for Western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat MRPS16. Suitable for mitochondrial translation studies and COXPD2 research.

[Learn More](#)

Parameter	Value
Product Name	MRPS16 Rabbit pAb
Remarks/Alias	S16mt; COXPD2; RPMS16; CGI-132; MRP-S16; MRPS16
Species	Human
GeneID	51021
Immunogen	Recombinant protein corresponding to amino acids 1-137 of human MRPS16 (NP_057149.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3D3
Protein Weight	15kDa
Shipping	Ice bag