

Myh10 Rabbit Polyclonal Antibody For Wb & Elisa - P35580

Item Number: CM0015151



Introduction

Rabbit polyclonal antibody against human MYH10 (Myosin-10), validated for Western blot and ELISA. Detects human, mouse, and rat MYH10 with a predicted molecular weight of 229 kDa. Ideal for studies on cell motility, cytokinesis, and viral entry.

[Learn More](#)

Parameter	Value
Product Name	MYH10 Rabbit pAb
Remarks/Alias	NMMHCB; NMMHC-IIB; MYH10
Species	Human
Gene ID (Human)	4628
Gene ID	4628
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1908-2007 of human MYH10 (NP_001242941.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35580
Protein Weight	229kDa
Shipping	Ice bag