

Ncapd2 Polyclonal Antibody For Wb And Elisa - Q15021

Item Number: CM0022552

Introduction

NCAPD2 rabbit polyclonal antibody validated for WB and ELISA applications, with reactivity across human, mouse, and rat samples. Condensin complex subunit 1, essential for chromosome condensation and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	NCAPD2 Rabbit pAb
Remarks/Alias	CNAP1; CAP-D2; MCPH21; hCAP-D2; NCAPD2
Species	Human
Gene ID (Human)	9918
Gene ID	9918
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1250-1350 of human NCAPD2 (NP_055680.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB,ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15021
Protein Weight	157kDa
Shipping	Ice bag