


Nfs1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y697

Item Number: CM0018225



Introduction

Rabbit polyclonal antibody targeting NFS1 cysteine desulfurase. Validated for WB, IF/ICC, and ELISA in human, mouse, and rat samples. Suitable for iron-sulfur cluster biogenesis and mitochondrial research.

[Learn More](#)

Parameter	Specification
Product Name	NFS1 Rabbit pAb
Remarks/Alias	IscS; NIFS; COXPD52; HUSSY-08; NFS1
Species	Human
GeneID (Human)	9054
GeneID	9054
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 208-457 of human NFS1 (NP_066923.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y697
Protein Weight	50kDa
Shipping	Ice bag