


Nr5A2 Rabbit Pab For Wb, If-P, Ihc-P, Elisa - O00482

Item Number: CM0027438



Introduction

Rabbit polyclonal antibody against human NR5A2 (LRH-1/FTF). Validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for research on cholesterol metabolism, embryonic development, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	NR5A2 Rabbit pAb
Remarks/Alias	B1F; CPF; FTF; B1F2; LRH1; LRH-1; FTZ-F1; hB1F-2; FTZ-F1beta; NR5A2
Species	Human
Gene ID (Human)	2494
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-541 of human NR5A2 (NP_995582.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00482
Protein Weight	61kDa
Shipping	Ice bag