

Orexin B Rabbit Monoclonal Antibody For Wb, Ihc, If/Icc, Elisa - O43612

Item Number: CM0007205



Introduction

Rabbit monoclonal antibody against human Orexin B/Hypocretin-2. Validated for WB, IHC-P, IF/ICC, ELISA. Ideal for studying neuropeptide regulation of feeding, sleep-wake cycles, and energy metabolism.

[Learn More](#)

Parameter	Value
Product Name	Orexin B Rabbit mAb
Remarks/Alias	OX; PPOX; NRCLP1; Orexin B
Species	Human
Gene ID (Human)	3060
Gene ID	3060
Immunogen	Recombinant protein of human Orexin B/Hypocretin-2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O43612
Protein Weight	13kDa
Shipping	Ice bag