


Pan Cadherin Rabbit Monoclonal Antibody (Mab) For Wb, Ihc-P, Elisa - P12830 / P19022 / P22223 / P33151 / P55283 / P39038

Item Number: CM0009306



Introduction

High-quality pan cadherin rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Recognizes human and mouse cadherins. Ideal for cell adhesion and epithelial biology research.

[Learn More](#)

Parameter	Specification
Product Name	Pan Cadherin Rabbit mAb
Remarks/Alias	Arc-1; CD324; CDHE; ECAD; LCAM; UVO; Pan Cadherin
Species	Human
Gene ID (Human)	999/1000/1001/1002
Gene ID	999/1000/1001/1001/1002/12561
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P12830/P19022/P22223/P33151/P55283/P39038
Protein Weight	140 kDa
Shipping	Ice bag