

Pan-Methyl-Histidine Rabbit Polyclonal Antibody For Wb, Db, Elisa Applications

Item Number: CM0000207

Introduction

This pan-methyl-histidine rabbit polyclonal antibody, raised against a synthetic methylated histidine peptide, is suitable for WB, DB, and ELISA assays. Detects pan-methyl-histidine modifications in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	pan-methyl-Histidine Rabbit pAb
Gene ID (Human)	pan Hme
Immunogen	Synthetic peptide of Methylated Histidine (His)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag