

Pan-Trimethyl-Lysine (Pan Kme3) Rabbit Polyclonal Antibody For Wb And Elisa

Item Number: CM0000210

Introduction

Rabbit polyclonal antibody targeting pan-trimethyl-lysine (pan Kme3) for WB and ELISA applications. Exhibits cross-reactivity with human, mouse, and rat samples. Suitable for epigenetic and post-translational modification research.

[Learn More](#)

Parameter	Value
Product Name	pan-TriMethyl-lysine Rabbit pAb
Gene ID (Human)	pan Kme3
Immunogen	Synthetic peptide containing trimethylated K
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag