

Pd-L1/Cd274 Rabbit Monoclonal Antibody For Fc, Elisa, Inhibition - Q9Nzq7

Item Number: CM0003712



Introduction

Rabbit monoclonal antibody against human PD-L1/CD274 (B7-H1), validated for flow cytometry, ELISA, and inhibition assays. Recombinant immunogen (aa 1-290), 33 kDa, for immune checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	PD-L1/CD274 Rabbit mAb
Remarks/Alias	B7-H, B7H1, PDL1, PD-L1, hPD-L1, PDCD1L1, PDCD1LG1, PD-L1/CD274
Species	Human
GeneID (Human)	29126
GeneID	29126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PD-L1/CD274 (Q9NZQ7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA, Inhibition
Cross Reactivity	Human
SWISS	Q9NZQ7
Protein Weight	33kDa
Shipping	Ice bag