

Perk Rabbit Polyclonal Antibody For If/Icc, Elisa - Q9Nzj5

Item Number: CM0024335



Introduction

Rabbit polyclonal antibody targeting human PERK/EIF2AK3, validated for IF/ICC and ELISA. Cross-reacts with mouse. Recognizes the ER stress sensor kinase involved in unfolded protein response and eIF2 α phosphorylation. Ideal for cancer, metabolic and stress research.

[Learn More](#)

Parameter	Value
Product Name	PERK Rabbit pAb
Remarks/Alias	PEK; WRS; PERK
Species	Human
GeneID (Human)	9451
GeneID	9451
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-600 of human PERK (NP_004827.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NZJ5
Protein Weight	125kDa
Shipping	Ice bag