

Phospho-Akt1-S129 Rabbit Pab For Wb, Elisa - P31749

Item Number: CM0000989



Introduction

Phospho-AKT1-S129 Rabbit pAb is a polyclonal antibody raised in rabbit, validated for WB and ELISA. It specifically recognizes phosphorylated serine 129 of human AKT1, suitable for studying AKT signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1-S129 Rabbit pAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-S129
Species	Human
GeneID (human)	207
GeneID	207
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S129 of human Akt1 (NP_005154.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag