


Phospho-C-Jun-T93 Rabbit Pab For Wb, Elisa - P05412

Item Number: CM0000487



Introduction

Phospho-specific rabbit polyclonal antibody targeting phosphorylated Thr93 of human c-Jun (AP-1). Validated for WB and ELISA; reacts with human and mouse. Ideal for studying AP-1 transcription factor activity and signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Jun-T93 Rabbit pAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-T93
Species	Human
GeneID (human)	3725
GeneID	3725
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T93 of human JUN (NP_002219.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag