

Phospho-C/Ebpb-T235 Rabbit Polyclonal Antibody For Wb, Elisa - P17676

Item Number: CM0000676



Introduction

Phospho-C/EBPB-T235 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and mouse C/EBPB phosphorylated at T235, a key transcription factor in immune response and adipogenesis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-C/EBPB-T235 Rabbit pAb
Remarks/Alias	TCF5; IL6DBP; NF-IL6; C/EBP-beta; Phospho-C/EBPB-T235
Species	Human
Gene ID (Human)	1051
Gene ID	1051
Immunogen	A synthetic phosphorylated peptide around T235 of human C/EBPB (NP_005185.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17676
Protein Weight	36kDa
Shipping	Ice bag