


Phospho-Chk1-S317 Rabbit Polyclonal Antibody (Wb, Elisa) - O14757

Item Number: CM0000628



Introduction

Highly specific rabbit polyclonal antibody recognizing endogenous levels of Chk1 phosphorylated at Ser317. Validated for Western blotting (WB) and ELISA in human samples. Ideal for cell cycle and DNA damage research.

[Learn More](#)

Parameter	Details
Product Name	Phospho-Chk1-S317 Rabbit pAb
Remarks/Alias	CHK1; Phospho-Chk1-S317
Species	Human
Gene ID (Human)	1111
Gene ID	1111
Immunogen	A synthetic phosphorylated peptide around S317 of human Chk1 (NP_001265.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag