


Phospho-Chk1-S345 Rabbit Mab For Wb, Elisa - O14757

Item Number: CM0000505



Introduction

Rabbit monoclonal antibody detecting endogenous phosphorylated Chk1 at Ser345, validated for western blot and ELISA with human and mouse cross-reactivity, targeting the key cell cycle checkpoint kinase.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Chk1-S345 Rabbit mAb
Remarks/Alias	CHK1; Phospho-Chk1-S345
Species	Human
GeneID (Human)	1111
GeneID	1111
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S345 of human Phospho-Chk1-S345 (NP_001107593.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag