


Phospho-Gab1 (Tyr659) Rabbit Polyclonal Antibody For Wb, Elisa - Q13480

Item Number: CM0001009



Introduction

Rabbit polyclonal antibody detecting phosphorylated GAB1 at tyrosine 659. Validated for Western blot (WB) and ELISA; reacts with human and mouse. Ideal for signaling pathway studies involving EGFR, INSR, and MET/HGF. Mol weight ~77 kDa.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-GAB1-Y659 Rabbit pAb
Remarks/Alias	DFNB26; Phospho-GAB1-Y659
Species	Human
Gene ID (Human)	2549
Gene ID	2549
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y659 of human GAB1 (NP_002030.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13480
Protein Weight	77kDa
Shipping	Ice bag