


Phospho-Grin2B-Y1474 Rabbit Polyclonal Antibody For Wb, Elisa - Q13224

Item Number: CM0001029



Introduction

Rabbit polyclonal antibody targeting phosphorylated GRIN2B at Y1474. Validated for WB and ELISA with human samples. Useful for studying NMDA receptor signaling, synaptic plasticity, and neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GRIN2B-Y1474 Rabbit pAb
Remarks/Alias	NR3, MRD6, NR2B, hNR3, DEE27, EIEE27, GluN2B, NMDAR2B, Phospho-GRIN2B-Y1474
Species	Human
GeneID (Human)	2904
GeneID	2904
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around Y1474 of human GRIN2B (NP_000825.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13224
Protein Weight	166kDa
Shipping	Ice bag