


Phospho-Histone H2Ax-S139 Rabbit Pab For Wb, If/Icc, Ihc-P, Elisa - P16104

Item Number: CM0000512



Introduction

Phospho-Histone H2AX-S139 Rabbit pAb detects human, mouse, rat targets. Designed for DNA damage, repair, and cell cycle studies via WB, IF/ICC, IHC-P, ELISA. Research-grade antibody.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H2AX-S139 Rabbit pAb
Remarks/Alias	H2A.X; H2A/X; H2AFX; Phospho-Histone H2AX-S139
Species	Human
GeneID (Human)	H2A
GeneID	3014
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S139 of human Histone H2AX (NP_002096.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16104
Protein Weight	15 kDa
Shipping	Ice bag