


Phospho-Histone H3-T32 Rabbit Pab For Wb, Ihc-P, Elisa - Q16695 / P68431

Item Number: CM0000649



Introduction

Phospho-Histone H3-T32 Rabbit pAb (AP0897) is a rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets phosphorylated Thr32 of histone H3. Ideal for chromatin and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3-T32 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Phospho-Histone H3-T32
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag