


Phospho-Hsf1-S326 Rabbit Pab For Wb, Elisa - Q00613

Item Number: CM0001873



Introduction

Rabbit polyclonal antibody targeting human HSF1 phosphorylated at serine 326. Suitable for western blot and ELISA applications in heat shock response, transcription regulation, and stress research. Detects a 57 kDa band.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSF1-S326 Rabbit pAb
Remarks/Alias	HSTF1; Phospho-HSF1-S326
Species	Human
GeneID	3297
Immunogen	A synthetic phosphorylated peptide around S326 of human HSF1 (NP_005517.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q00613
Protein Weight	57kDa
Shipping	Ice bag