

Phospho-Insr-Y1146/Igf1R-Y1131 Rabbit Pab For Wb, Elisa - P06213 / P08069

Item Number: CM0001126



Introduction

Rabbit polyclonal antibody specifically recognizes phospho-INSR (Y1146) and phospho-IGF1R (Y1131) in human samples. Validated for Western blot and ELISA, it is an essential tool for insulin/IGF signaling research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-INSR-Y1146/IGF1R-Y1131 Rabbit pAb
Remarks/Alias	INSR/IGF1R; Phospho-INSR-Y1146/IGF1R-Y1131
Species	Human
GeneID (Human)	3480/3643
GeneID	3643/3480
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around Y1146 of human IGF1RINSR (NP_000199.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06213/P08069
Protein Weight	155kDa
Shipping	Ice bag