

Phospho-Ire1-S724 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75460

Item Number: CM0001061



Introduction

Phospho-specific rabbit monoclonal antibody targeting IRE1 (ERN1) phosphorylated at Ser724. Validated for WB, IHC-P, IF/ICC, and ELISA. Reliable tool for investigating ER stress and UPR pathways in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRE1-S724 Rabbit mAb
Remarks/Alias	IRE1; IRE1P; IRE1a; hIRE1p; Phospho-IRE1-S724
Species	Human, Mouse, Rat
Gene ID (Human)	2081
Gene ID	2081
Immunogen	Synthetic phosphorylated peptide around S724 of human IRE1 (NP_001424.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75460
Protein Weight	110kDa
Shipping	Ice bag