

Phospho-Map2K1-S217/Map2K2-S221 Rabbit Polyclonal Antibody For Wb, Elisa - Q02750 / P36507

Item Number: CM0001753



Introduction

Phospho-specific rabbit polyclonal antibody detecting MAP2K1 (phospho S217) and MAP2K2 (phospho S221). Suitable for WB and ELISA in human, mouse, and rat samples. Ideal for studying MAPK/ERK signaling cascade activation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP2K1-S217/MAP2K2-S221 Rabbit pAb
Remarks/Alias	MEK1/MEK2; Phospho-MAP2K1-S217/MAP2K2-S221
Species	Human
Gene ID (Human)	5604/5605
Immunogen	A synthetic phosphorylated peptide around S217 & S221 of human MAP2K1/MAP2K2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750/P36507
Protein Weight	40kDa/43kDa/44kDa
Shipping	Ice bag