

Phospho-Marcks-S162 Rabbit Polyclonal Antibody For Wb, Elisa - P29966

Item Number: CM0001025

Introduction

Rabbit polyclonal antibody specific for phosphorylated MARCKS (Ser162). Validated for WB and ELISA, reacts with human, mouse, and rat. Ideal for studying PKC signaling, cytoskeletal dynamics, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MARCKS-S162 Rabbit pAb
Remarks/Alias	MACS; 80K-L; PKCSL; PRKCSL; Phospho-MARCKS-S162
Species	Human
GeneID (Human)	4082
GeneID	4082
Immunogen	Synthetic peptide A phospho specific peptide corresponding to residues surrounding S162 of human MARCKS
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29966
Protein Weight	32kDa
Shipping	Ice bag