

Phospho-Mdm2-S166 Rabbit Mab For Ihc-P And Elisa - Q00987

Item Number: CM0001665



Introduction

Phospho-MDM2-S166 Rabbit mAb specifically detects human MDM2 phosphorylated at serine 166. Validated in IHC-P and ELISA; cross-reacts with mouse and rat. Supports research on the p53-MDM2 axis and apoptosis regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MDM2-S166 Rabbit mAb
Remarks/Alias	HDMX; LSKB; hdm2; ACTFS; Phospho-MDM2-S166
Species	Human
Gene ID (Human)	4193
Gene ID	4193
Immunogen	Synthetic peptide: A phospho synthetic peptide corresponding to residues surrounding S166 of Human MDM2.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00987
Protein Weight	55kDa
Shipping	Ice bag