

Phospho-Met-Y1349 Rabbit Polyclonal Antibody For Wb, Elisa - P08581

Item Number: CM0000548



Introduction

Validated Phospho-MET (Y1349) rabbit polyclonal antibody for WB and ELISA, cross-reactive with human and rat. Detects hepatocyte growth factor receptor phosphorylation, essential for studying MET signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MET-Y1349 Rabbit pAb
Remarks/Alias	DA11; HGFR; AUTS9; RCCP2; c-Met; DFNB97; Phospho-MET-Y1349
Species	Human
GeneID (Human)	4233
GeneID	4233
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y1349 of human MET (NP_000236.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P08581
Protein Weight	156kDa
Shipping	Ice bag