


Phospho-Rb-S795 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P06400

Item Number: CM0001527



Introduction

Phospho-RB-S795 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Recognizes human, mouse, and rat RB1 phosphorylated at Ser795. Essential tool for cell cycle and tumor suppressor research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-RB-S795 Rabbit mAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1; Phospho-RB-S795
Species	Human
Gene ID	5925
Immunogen	A synthetic phosphorylated peptide around S795 of human RB.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06400
Protein Weight	106kDa
Shipping	Ice bag