


Phospho-S6 Ribosomal Protein (Rps6)-S235 Rabbit Polyclonal Antibody For Wb, Elisa - P62753

Item Number: CM0001308



Introduction

Phospho-specific rabbit polyclonal antibody targeting RPS6 phosphorylated at Ser235. Validated for WB and ELISA, with cross-reactivity in human, mouse, and rat. Ideal for studying cell growth and proliferation signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-S6 Ribosomal Protein (RPS6)-S235 Rabbit pAb
Remarks/Alias	S6; eS6; Phospho-S6 Ribosomal Protein (RPS6)-S235
Species	Human
Gene ID (Human)	6194
Immunogen	A phospho specific peptide corresponding to residues surrounding S235 of human RPS6
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62753
Protein Weight	29kDa
Shipping	Ice bag