

## Phospho-Stat2-Y690 Rabbit Polyclonal Antibody For Wb And Elisa - P52630

Item Number: CM0001481



### Introduction

Rabbit polyclonal antibody recognizing phosphorylated STAT2 at Tyr690. Suitable for Western blot (WB) and ELISA. Reacts with human, mouse, and rat. Ideal for studying interferon signaling, antiviral responses, and mitochondrial dynamics.

[Learn More](#)

| Parameter        | Value                                                       |
|------------------|-------------------------------------------------------------|
| Product Name     | Phospho-STAT2-Y690 Rabbit pAb                               |
| Remarks/Alias    | P113; IMD44; ISGF-3; PTORCH3; STAT113; Phospho-STAT2-Y690   |
| Species          | Human                                                       |
| Gene ID (Human)  | 6773                                                        |
| Gene ID          | 6773                                                        |
| Immunogen        | Synthetic phosphorylated peptide around Y690 of human STAT2 |
| Source           | Rabbit                                                      |
| Category         | Phosphorylated Antibodies                                   |
| Application      | WB, ELISA                                                   |
| Cross Reactivity | Human, Mouse, Rat                                           |
| SWISS            | P52630                                                      |
| Protein Weight   | 98kDa                                                       |
| Shipping         | Ice bag                                                     |