


Phospho-Stk11-S334 Rabbit Polyclonal Antibody For Wb, Elisa - Q15831

Item Number: CM0001431



Introduction

Phospho-STK11-S334 Rabbit Polyclonal Antibody detects human STK11 phosphorylated at Ser334, with cross-reactivity to mouse and rat. Suitable for WB and ELISA; ideal for studying tumor suppression, cell metabolism, and DNA damage.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STK11-S334 Rabbit pAb
Remarks/Alias	PJS; LKB1; hLKB1; Phospho-STK11-S334
Species	Human
GeneID	6794
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S334 of human STK11 (NP_000446.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15831
Protein Weight	49kDa
Shipping	Ice bag