

Phospho-Tau-T217 Rabbit Mab For Wb, Elisa - P10636-8

Item Number: CM0000670

Introduction

Rabbit recombinant monoclonal antibody against phospho-Tau (Thr217) validated for Western blot and ELISA. Specifically detects human tau, aiding Alzheimer's disease and tauopathy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T217 Rabbit mAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T217
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	A synthetic phosphorylated peptide around T217 of human Tau
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10636-8
Protein Weight	33-46kDa/79-81kDa
Shipping	Ice bag