


Pias1 Rabbit Polyclonal Antibody For Wb And Elisa - 075925

Item Number: CM0013062



Introduction

PIAS1 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, and rat. Targets E3 SUMO-protein ligase PIAS1 (075925), a key transcriptional coregulator in the STAT, p53, and hormone signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PIAS1 Rabbit pAb
Remarks/Alias	GBP; ZMIZ3; DDXBP1; GU/RH-II; PIAS1
Species	Human
GeneID (human)	8554
GeneID	8554
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human PIAS1 (NP_057250.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075925
Protein Weight	72kDa
Shipping	Ice bag