

Ppl Rabbit Polyclonal Antibody For Wb, If-P, Elisa - O60437

Item Number: CM0018080



Introduction

Perioplakin (PPL) rabbit polyclonal antibody validated for WB, IF-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studying cornified envelope and desmosome biology.

[Learn More](#)

Parameter	Value
Product Name	PPL Rabbit pAb
Remarks/Alias	PPL
Species	Human
GeneID (Human)	5493
GeneID	5493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human PPL (NP_002696.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60437
Protein Weight	205kDa
Shipping	Ice bag