


Prkar1A Rabbit Polyclonal Antibody (Ko Validated) For Wb, If/Icc, Elisa - P10644

Item Number: CM0016684



Introduction

High-quality rabbit polyclonal antibody against human PRKAR1A, validated by knockout, suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Research use for cAMP signaling pathway.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PRKAR1A Rabbit pAb
Remarks/Alias	CAR; CNC; CNC1; PKR1; TSE1; ADOHR; PPNAD1; PRKAR1; ACRDYS1; 1A
Species	Human
GeneID (Human)	5573
GeneID	5573
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human PRKAR1A (NP_002725.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10644
Protein Weight	43kDa
Shipping	Ice bag