

Prkra/Pact Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - 075569

Item Number: CM0004133



Introduction

Recombinant rabbit monoclonal antibody targeting human PRKRA/PACT (075569). Validated for Western blot, IHC-P, and ELISA; cross-reacts with mouse and rat. Suitable for studies on PKR signaling, RNA silencing, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	PRKRA/PACT Rabbit mAb
Remarks/Alias	RAX; PACT; DYT16; HSD14; PRKRA/PACT
Species	Human
GeneID (Human)	8575
GeneID	8575
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 214-313 of human PRKRA/PACT (075569).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075569
Protein Weight	34kDa
Shipping	Ice bag