

Pseudouridine (5-Ribosyluracil) Rabbit Polyclonal Antibody For Dot Blot And Elisa

Item Number: CM0000042

Introduction

Rabbit polyclonal antibody against pseudouridine/5-ribosyluracil, a modified nucleoside. Validated for dot blot and ELISA with species-independent reactivity. Ideal for epitranscriptomic and diagnostic research applications.

[Learn More](#)

Parameter	Value
Product Name	Pseudouridine / 5-ribosyluracil Rabbit pAb
Immunogen	Chemical compounds corresponding to Pseudouridine / 5-ribosyluracil / Y
Source	Rabbit
Category	Other Modified Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag