

Rabbit Anti-Bgat (Abo) Polyclonal Antibody For Wb And Elisa - P16442

Item Number: CM0023913



Introduction

Rabbit polyclonal antibody against human BGAT (ABO transferase) validated for WB and ELISA. Reacts with human, mouse, and rat samples. Ideal for ABO blood group antigen research.

[Learn More](#)

Parameter	Value
Product Name	BGAT Rabbit pAb
Remarks/Alias	GTB; NAGAT; A3GALNT; A3GALT1; BGAT
Species	Human
GeneID	28
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 54-354 of human BGAT (NP_065202.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16442
Protein Weight	41kDa
Shipping	Ice bag